

Medonic M-series Blood Cell Analyzer



Open Vial Sampling System

Cap Piercing Sampling System

Patent Finger Stick Technology

20ul draw is all it takes. Eliminates the need for costly drawing devices

Automate Your Records & Data

Integrate a Data Management Package with your system.

- Ensures compliance with all patient regulations regarding patient data, patient records, quality control and general record keeping.
- Easy to use
- Affordable solution for the physician`s office.
- Easy interface to EMR

Specifications

16-Parameter, 3 Part Diff
Up to 60 Samples Per Hour
Results in 45 Seconds
Micro Sampling
LCD Color Touch Screen Display
Histogram on Screen Display
On Board Step By Step Operation Instructions
Bar Code Reader
Sample Memory
Control Memory
Program Files
Clot Removal
On Board QC Including LJ, XB
Blood Detectors
Auto Hgb Blank
Hgb Correction
Reagent Level Sensors
Warning Flags
Parameter Flags
Interface Ready
Zero Daily Maintenance

What Differentiates The M Series

Accurate and Precise - Delivers excellent patient results.

Bar Coding - Standard with all systems giving you positive patient identification and reduces transcription errors.

Zero Daily Maintenance - No daily shutdown, engineered to eliminate lengthy cleaning procedures.

Service - Highly reliable, excellent uptime record.

Options

Key Board
Data Management Package

Small Footprint

16.15 x 11.43 x 18.12

Weight

39.6 Pounds - Open / Closed Vial
Power Consumption 100VA
Power Consumption (stand-by) 20VA
Frequency 50/60Hz